

BNP Monitoring


Fingerstick POC Platform – BNP Test System



- Easy to use
- Connectivity enabled – cellular transmission to data center
- Ability to capture other data points


Device View – Stirling BNP Test Strip

- Low Volume Fingertstick
- 12uL total blood volume
- ~15 minute TTR
- On Board Controls
- Electrochemical Immunoassay
- Biomarker Platform
- Potential to add further content (multi-markers) in next generation




System Views – Stirling BNP Test System

Top View



Rear View



ROM key port
(ROM key not
inserted)


Front View



RS232 data
port connector Mains power
connector

Back View

System Views – PST



System Views – PST

